
# MASSERANN

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
## INTRODUCTION

Rotary NITI systems are at the origin of a new demand for devices to remove broken instruments. The Masserann Micro Kit will be the safety parachute in your dental practice.

## PROTOCOL

### 1/ Open access

Enlarge canal orifice using a round bur.

Use Gates drill to create straight access to the fragment.

### 2/ Clear fragment

Before use, make sure that the trepan is in good condition and verify the quality of the tines' cutting edges.

Use trepan bur to set free approximately 4 mm of fragment (slow counter-clockwise rotation). Masserann trepan burs can be used either by hand with the short manual handle or with a reducing contra-angle (300-600 rpm) with attachment system according to ISO 1797-1 standard (tightly screw with the flat spanner).

### 3/ Grip and remove fragment

Slide the opened extractor over the freed end of the fragment.

Hold the extractor firmly in place and turn the screw head until the fragment is gripped.

Remove the fragment:

- In counter-clockwise rotation for all cutting root canal instruments.
- In clockwise rotation for all filling instruments (ex: spiral fillers).

Your Endo Specialist™

**Bear in mind that....**



Fragments must always be removed in rotation – never pull out only.  
 Direct straight access is a must: Masserann tools are straight.  
 Positioning the instruments is most important:  
 Use your finger to guide the trepan in the right direction.  
 Two X-rays from different angles will help checking the orientation of the fragment before and when applying the extractor.

**Cleaning:**

Immerge in a detergent bath, possibly ultrasonic cleaning bath. Cleaning, disinfection and sterilization instructions for root canal instruments recommended by MICRO-MEGA® also apply to the Masserann instruments.  
 Sterilization: autoclave - 135°C/2.1 bars (275°F/29psi) – 18 min minimum under sterilization bag. For sterilization, we recommend that the tray is placed inside the lid (the tray is open). The tray must be sterilized with the instruments inside.

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